



# FreeStyle

## Optium Neo

Blood Glucose and Ketone Monitoring System

### Owner's Setup Guide User's Manual




#### TABLE OF CONTENTS

- 1 FreeStyle Optium Neo At-A-Glance**
- 2 Meter Symbols**
- 3 Quick Start - Testing Your Blood Glucose**
- 4 Intended Use**
- 5 Getting to Know Your Meter**  
Turning Your Meter On and Off  
Checking Your Meter Screen Every Time You Turn It On
- 6 Setting Up the Meter**
- 7 Testing Your Blood Glucose or Blood Ketone**  
Prepare to Test  
Performing a Blood Glucose or Blood Ketone Test
- 8 Understanding Blood Glucose Test Results and Patterns**  
Low Blood Glucose Results  
High Blood Glucose Results
- 9 Understanding Blood Ketone Test Results**
- 10 View Logbook**  
Viewing Your Logbook Events  
Viewing Your Blood Glucose Averages

#### TABLE OF CONTENTS

- 11 Glucose and Ketone Control Solutions**  
Performing a Control Solution Test
- 12 Transferring Meter Data to a Computer**
- 13 Error Messages**
- 14 Taking Care of Your Meter**  
Changing Batteries  
Cleaning Your Meter
- 15 Troubleshooting**
- 16 Insulin Dose Logging**  
Set Up for Insulin Dose Logging  
How to Use Insulin Dose Logging  
How to Log Additional Rapid-Acting Insulin Doses
- 17 Meter Specifications**
- 18 Other Symbols**
- 19 References**


#### 1 FreeStyle Optium Neo At-A-Glance



#### 2 Meter Symbols

Symbol	What It Means	Symbol	What It Means
	Logbook		Breakfast insulin dose
	Meter ready for sample application		Lunch insulin dose
	Control solution result		Dinner insulin dose
	Units of insulin		Fasting test
	Insulin	<b>KET</b>	Ketone
	Morning long-acting insulin dose		Low battery
	Evening long-acting insulin dose		Connected to computer
	Insulin dose setup		Setup mode
	Low blood glucose or pattern		High blood glucose or pattern

#### 3 Quick Start – Testing Your Blood Glucose



#### 4 Intended Use

The FreeStyle Optium Neo Blood Glucose and Ketone Monitoring System is for use outside the body only (*in vitro* diagnostic use) for self testing or professional use as an aid in the management of diabetes.

Use for measuring glucose in fresh whole blood samples taken from fingers, forearm, upper arm or the base of the thumb. Use for measuring ketone ( $\beta$ -hydroxybutyrate) in fresh whole blood samples from fingers only.

The system may also aid in the management of diabetes by providing the user with suggested recommendations to insulin dose(s) based on healthcare professional entered data.

##### IMPORTANT:

- Use only FreeStyle Optium blood glucose test strips and FreeStyle Optium blood  $\beta$ -Ketone test strips. Other test strips may produce inaccurate results.
- See test strip instructions for use for more information about sample types.
- Read the instructions in this Owner's Setup Guide. Failure to follow instructions may cause incorrect results. Practice the testing procedures before using the meter.
- Follow your healthcare professional's advice when testing blood glucose levels and blood ketone levels.
- Observe caution when using around children. Small parts may constitute a choking hazard.

#### 5 Getting to Know Your Meter

##### Turning Your Meter On and Off

To turn your meter on:


- Press the , or
- Insert a strip

To turn your meter off:

- Press and hold the for 3 seconds, or
- Do nothing for 2 minutes

##### Checking Your Meter Screen Every Time You Turn It On

###### Meter Off Screen



Start-up Test Screen

Your meter screen should be fully black when powered off. Each time you turn on your meter, a white start-up test screen will appear for 1 second.



If you see any *white* segments in the black off screen, or any *black* segments in the white test screen, there may be a problem with the meter. Contact Customer Care.

*Note: If the meter battery is low, the will appear in both the meter off screen and start-up test screen.*


#### 6 Setting Up the Meter

To confirm that the date and time are set correctly, follow the steps below.

##### Set Time




##### Set Date



#### 10 View Logbook

Your meter logbook can store up to 1,000 events – including blood glucose, ketone, and control solution results, insulin doses, and other meter information.

##### Viewing Your Logbook Events



#### 7 Testing Your Blood Glucose or Blood Ketone

##### IMPORTANT:

- Only use a test strip once.
- Read the test strip instructions for use before performing your first blood glucose or ketone test. It contains important information and will tell you how to store and handle the test strips.
- The meter and its accessories are for use by a single person.
- Refer to the lancing device insert for detailed instructions on how to use the lancing device.
- Do not put urine on the test strip.

##### Check your ketone:

- When you have an illness
- When your blood glucose is above 13.3 mmol/L
- When you and your healthcare professional feel it is necessary

##### Prepare to Test

###### 1. Select your test site.

Test site choices for blood glucose testing are fingers, forearm, upper arm, or base of the thumb. Use only fingertip blood samples for blood ketone testing.

*Note: Avoid moles, veins, bones, and tendons. Bruising may occur at the test site. If bruising occurs, consider selecting another site.*

##### Do not use blood samples from alternative sites when:

- You think your blood glucose is low or changing rapidly
- You have been diagnosed with hypoglycaemic unawareness
- Alternative site results do not match the way you feel
- You are within two hours of eating a meal, taking insulin, or exercising




###### 2. Wash your hands and the test site with soap and warm water.

- Rinse and dry thoroughly.
- Do not use lotion or cream on the test site.

###### 3. Check test strip expiry date.

Do not use expired test strips; they may cause inaccurate results.

##### Performing a Blood Glucose or Blood Ketone Test



#### 8 Understanding Blood Glucose Test Results and Patterns

The meter displays blood glucose results in mmol/L. The unit of measurement is preset. You cannot change this setting.

**IMPORTANT:** The meter displays results from 1.1 – 27.8 mmol/L. Low or high blood glucose results can indicate a potentially serious medical condition.

Low Blood Glucose Results	High Blood Glucose Results
<p><b>If You See ...</b></p> <p></p> <p>A solid red arrow</p> <p></p> <p>A blinking red arrow</p> <p><b>LO</b></p> <p>Appears when result is lower than 1.1 mmol/L</p>	<p><b>If You See ...</b></p> <p></p> <p>A solid yellow arrow</p> <p></p> <p>A blinking yellow arrow</p> <p><b>HI</b></p> <p>Appears when result is higher than 27.8 mmol/L</p>
<p><b>What It Means</b></p> <p>Appears when result is lower than 3.9 mmol/L or the target set on the meter by your healthcare professional.</p> <p>A pattern of low glucose has developed. If 2 low results occur within the past 5 days AND both are within the same 3-hour time period, the meter will display a blinking .</p> <p>Severe low blood glucose or There may be a problem with the test strip.</p>	<p><b>What It Means</b></p> <p>Appears when result is higher than 13.3 mmol/L or the target set on the meter by your healthcare professional.</p> <p>A pattern of high glucose has developed. If 3 high results occur within the past 5 days AND all are within the same 3-hour time period, the meter will display a blinking .</p> <p>Blood glucose level is higher than or equal to 13.3 mmol/L.</p> <p>Severe high blood glucose or There may be a problem with the test strip.</p>
<p><b>What To Do</b></p> <p>Follow your healthcare professional's advice to treat low blood glucose.</p> <p>Follow your healthcare professional's advice to treat low blood glucose.</p> <p>Repeat the test with a new test strip. If the result is LO, contact your healthcare professional <b>immediately</b>.</p>	<p><b>What To Do</b></p> <p>Follow your healthcare professional's advice to treat high blood glucose.</p> <p>Follow your healthcare professional's advice to treat high blood glucose.</p> <p>Check blood ketone if checking ketones is part of your diabetes management program.</p> <p>Repeat the test with a new test strip. If the result is HI, contact your healthcare professional <b>immediately</b>.</p>

*Note: If you see the error messages E-3 or E-4, consult the Error Messages section in this Owner's Setup Guide.*

**IMPORTANT:** Contact your healthcare professional if you have symptoms that do not match your test result, and you have followed the instructions in this Owner's Setup Guide.

#### 9 Understanding Blood Ketone Test Results

The meter displays ketone results in mmol/L, from 0.0 – 8.0 mmol/L. The unit of measurement is preset. You cannot change this setting.

**IMPORTANT:** Follow your healthcare professional's advice before you make any changes to your diabetes management program.


Blood ketone is expected to be below 0.6 mmol/L.<sup>2</sup> High blood ketone may be caused by illness, fasting, vigorous exercise, or uncontrolled blood glucose levels.<sup>1,3</sup>

Repeat a blood ketone test using a new blood ketone test strip when:

- HI appears on the display
- Your result is unusually high
- You question your result
- Your blood ketone result is 0.0 mmol/L, but your blood glucose is higher than 16.7 mmol/L

Displays	What It Means	What To Do
Result is between 0.6 – 1.5 mmol/L	High blood ketone. A problem requiring medical assistance may be occurring.	Follow your healthcare professional's advice.
Result is higher than 1.5 mmol/L	You may be at risk of developing diabetic ketoacidosis (DKA). <sup>2,4</sup>	Contact your healthcare professional <b>immediately</b> .
<b>HI</b>	Very high blood ketone or There may be a problem with the test strip.	Repeat the test with a new test strip. If the result is HI, contact your healthcare professional <b>immediately</b> .

##### Viewing Your Blood Glucose Averages





## 11 Glucose and Ketone Control Solutions



A control solution test should be performed when you are not sure of your results and want to confirm that your meter and test strips are working properly.


### IMPORTANT:

- Use only MediSense glucose and ketone control solutions with the meter.
- Control solution results should fall within the control solution range printed on the test strip instructions for use.
- Check that the LOT number printed on the test strip foil packet and instructions for use match.
- Do not use control solution past its expiry date. Discard control solution 3 months after opening or on the expiry date printed on the bottle, whichever comes first. (Example: open April 15, discard July 15; write the discard date on the side of the bottle.)
- The control solution range is a target range for control solution only, not for blood glucose levels.
- Replace the cap securely on the bottle immediately after use.

### IMPORTANT: (continued)

- Do not add water or other liquid to control solution.
- Control solution results do not reflect your blood glucose.
- Contact Customer Care for information on how to obtain control solution.

### Performing a Control Solution Test



1. Open the foil test strip packet at the notch and tear down to remove the test strip.

2. Insert the test strip until the meter turns on.

#### Notes:

- Be sure to check that your meter screen is working properly each time you turn your meter on. If you see any white segments in the black off screen, or any black segments in the white test screen, there may be a problem with the meter. (See Section 5, Getting to Know Your Meter, for more details.)
- The meter turns off after 3 minutes of inactivity. Remove and reinsert the unused test strip to restart the meter.




The blinks, indicating the meter is ready for you to apply a sample to the test strip.

**Note:** KET will appear on the screen if you have inserted a purple blood ketone strip.



3. Press and hold the down arrow for 3 seconds to mark the test as a control solution test. The appears. The meter is now ready for you to apply control solution to the test strip.

**IMPORTANT:** The test result will be saved to memory as a blood result if not marked as a control solution test. This may affect your blood glucose averages.



### 4. Apply control solution to the test strip.

Shake the control solution bottle to mix the solution. Apply a drop of control solution to the white area at the end of the test strip in the area shown. The control solution is drawn into the test strip.



### 5. Hold the control solution to the test strip until:

- You see 3 short lines on the meter screen. This means you have applied enough control solution and the meter is reading the control solution.

#### Notes:

- If you are testing with a blood glucose strip, you will see a 5-second countdown. If you are testing with a blood ketone strip, you will see a 10-second countdown.
- Do not remove the test strip from the meter during the countdown.
- If the countdown does not start, remove and discard the used test strip, turn off the meter and try again with a new strip.

Examples:



### 6. View the result.

The test is complete (examples shown) when the result appears on the meter screen. The result is stored in memory as a control solution result.

Compare the control solution result to the range printed on the blood glucose or blood ketone test strip instructions for use. The result should fall within the range.

**Note:** KET appears with the result if performing a ketone control solution test.

#### Out of Range Control Solution Results:

- Repeat the test if control solution results are outside the range printed on the test strip instructions for use.
- Stop using the meter if control solution results are consistently outside the range printed on the test strip instructions for use. Contact Customer Care.

## 12 Transferring Meter Data to a Computer

Transferring meter data to a computer requires a compatible data management system. You will also need a micro USB cable to connect your meter's USB port to your computer. For more information, please contact Customer Care.

**WARNING:** To avoid the possibility of electric shock, never perform a blood glucose test while the meter is connected to the computer.

## 13 Error Messages

Message	What It Means	What To Do
E-1	The temperature is too hot or too cold for the meter to work properly	1. Move the meter and test strips to a location where the temperature is within the test strip operating range. (See test strip instructions for use for the appropriate range.) 2. Wait for the meter and test strips to adjust to the new temperature. 3. Repeat the test using a new test strip. 4. If the error reappears, contact Customer Care.
E-2	Meter error	1. Turn off the meter. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.
E-3	Blood drop is too small or Incorrect test procedure or There may be a problem with the test strip	1. Review the testing instructions. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.
E-4	The blood glucose level may be too high to be read by the system or There may be a problem with the test strip	1. Repeat the test using a new test strip. 2. If the error reappears, contact Customer Care.
E-5	Blood was applied to the test strip too soon	1. Review the testing instructions. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.

Message	What It Means	What To Do
E-6	Meter error	1. Check that you are using the correct strip for this meter. (See test strip instructions for use to verify your strip is compatible with this meter.) 2. Repeat the test using a test strip for use with your meter. 3. If the error reappears, contact Customer Care.
E-7	No coding required or Test strip may be damaged, used, or the meter does not recognise it	1. Check that you are using the correct test strip for this meter. (See test strip instructions for use to verify your strip is compatible with this meter.) 2. Repeat the test using a test strip for use with your meter. 3. If the error reappears, contact Customer Care.
E-9	Meter error	1. Turn off the meter. 2. Repeat the test using a new test strip. 3. If the error reappears, contact Customer Care.

## 14 Taking Care of Your Meter

### Changing Batteries

You will see this screen when your batteries are low.

**Note:** Your meter settings and logbook information will be saved when you change the batteries.

**IMPORTANT:** After you first see this warning, you can perform approximately 28 tests before you need to change the batteries.

**WARNING:** Batteries should be kept away from small children. If swallowed, contact a healthcare professional immediately.

Step	Action
1. Turn meter over and slide open the battery door on the side as shown.	

Step	Action
2. Remove the old batteries.	
3. Install new batteries with (+) facing up.	<b>Note:</b> The meter uses 2 replaceable CR 2032 coin cell batteries.
4. Slide the door back into place until it clicks.	<b>Note:</b> The next time you turn your meter on, it may prompt you to reset the time and date. (Refer to Section 6, Setting Up the Meter).

**Note:** When you no longer need the meter, remove the batteries and dispose of the batteries and meter in accordance with your local government regulations.

The European Battery Directive requires separate collection of spent batteries, aiming to facilitate recycling and to protect the environment. The batteries in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.

### Cleaning Your Meter

Step	Action
1	Clean the outside of the meter using a damp cloth and: <ul style="list-style-type: none"> <li>• Mild detergent/soap and water, or</li> <li>• 70% isopropyl alcohol, or</li> <li>• A mixture of 1 part household bleach, 9 parts water</li> </ul>
2	Allow the meter to dry.

**IMPORTANT:** Do not place the meter in water or other liquids. Avoid getting dust, dirt, blood, control solution, water, or any other substance in the meter's test strip port and USB port.

## 15 Troubleshooting

Problem	What It Means	What To Do
1. Test strip is inserted in the strip port and nothing happens.	Test strip is not inserted properly or fully into the meter	1. With the contact bars (3 black lines) facing up, insert the test strip into the meter until it stops. This turns on the meter. 2. If the meter still does not turn on, contact Customer Care.
	No batteries are installed; Batteries installed incorrectly	Refer to Section 14, Taking Care of Your Meter, on how to properly install batteries.
	Dead batteries	Change batteries. Reset date and time, if needed.
	Meter may be plugged into a computer (PC appears on meter screen)	Unplug the meter from the computer.
	Problem with the test strip	Try a new test strip.
	Problem with the meter	Contact Customer Care.
2. The test does not start after applying the blood sample.	Blood sample is too small	1. See test strip instructions for use for re-application instructions. 2. Repeat the test using a new test strip. 3. If the test still does not start, contact Customer Care.
	Sample applied after meter turns off	1. Review the testing instructions. 2. Repeat the test using a new test strip. 3. If the test still does not start, contact Customer Care.
	Problem with meter or test strip	1. Repeat the test using a new test strip. 2. If the test still does not start, contact Customer Care.

## 16 Insulin Dose Logging

**Introduction**  
This feature allows you to log insulin doses so they are recorded in the logbook. You can enable this feature at any time.

### Set Up for Insulin Dose Logging

1. While on the Home screen, press and hold for 3 seconds until the appears. Insulin Dose Logging is now enabled.

**Note:** To turn off this feature, repeat this step.

### How to Use Insulin Dose Logging

For long-acting insulin and/or insulin for breakfast, lunch, or dinner.

1. While on the Home screen, press .

2. Choose type of insulin dose.  
Using the chart below, press the button for the dose you want to log.

	Long-acting insulin
	Long-acting insulin
	Meal insulin
	Meal insulin
	Meal insulin

3. Press or to enter the actual dose amount.

4. Press to log dose.

means you took dose

### How to Log Additional Rapid-Acting Insulin Doses

(e.g. snacks, bedtime correction, etc.)

1. While on the Home screen, press and hold for 3 seconds until the screen changes.

means you took dose

2. Press or to enter the actual dose amount.

3. Press to log dose.

means you took dose

## 17 Meter Specifications

Assay method	Amperometry	Power source	Two CR 2032 lithium (coin cell) batteries
Automatic shutoff	At least two minutes of inactivity	Size	5.97 cm (w) x 8.68 cm (l) x 0.87 cm (d) 2.35 in (w) x 3.42 in (l) x 0.34 in (d)
Battery life	Up to 3000 tests	Storage temperature	-20 °C to 60 °C (-4 °F to 140 °F)
Measurement range	For blood glucose testing 1.1 - 27.8 mmol/L For blood ketone testing 0.0 - 8.0 mmol/L	System altitude	See test strip instructions for use
Memory	Up to 1000 events, including blood glucose, blood ketone, and control solution results, insulin doses, and other meter information	Weight	33 to 37 g (1.2 to 1.3 oz.) including batteries
Minimum computer requirements	System must only be used with EN60950-1 rated computers. Use a USB certified cable	<b>Note:</b> For test strip specifications, see test strip instructions for use.	
Operating relative humidity	10% to 90% (non-condensing)	Electromagnetic Compatibility (EMC): FreeStyle Optium Neo meter has been tested for both electrostatic discharge and radio frequency interference. Emissions are low and unlikely to interfere with other nearby electronic equipment. To limit radio frequency interference do not use the FreeStyle Optium Neo meter near cellular or cordless telephones, radio transmitters or other electrical or electronic equipment that are sources of electromagnetic radiation, as these may interfere with the proper operation of the meter. Avoid use of the device in very dry environments, as electrostatic discharges from synthetic materials (e.g., carpets) could cause damage.	
Operating temperature	Meter: 10 °C to 50 °C (50 °F to 122 °F) System: See test strip instructions for use		

## 18 Other Symbols

Symbol	What It Means	Symbol	What It Means
	Consult instructions for use		Caution
	Temperature limit		Use-by date
	Manufacturer		In vitro diagnostic medical device
	Batch code		Catalogue number
	Do not re-use		Recycle
	Date of manufacture		Serial number
	Do not drink		Sterilised using irradiation (lancets only)
	CE Mark		

The European Battery Directive requires separate collection of spent batteries, aiming to facilitate recycling and to protect the environment. The batteries in this product should be removed and disposed in accordance with local regulations for separate collection of spent batteries.

## 19 References

1. Schade DS, Eaton RP. Metabolic and clinical significance of ketosis. Special Topics in Endocrinology and Metabolism 1982; 4:1-27.
2. Wiggam MI, O'Kane MJ, Harper R, Atkinson AB, Hadden DR, Trimble ER, Bell PM. Treatment of diabetic ketoacidosis using normalization of blood 3-hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997; 20:1347-1352.
3. Harano Y, Kosugi K, Hyosu T, Suzuki M, Hidaka H, Kashiwagi A, Uno S, Shigeta Y. Ketone bodies as markers for Type 1 (insulin-dependent) diabetes and their value in the monitoring of diabetes control. Diabetologia 1984; 26:343-348.
4. Ukkata E. Diurnal variation of blood ketone bodies in insulin-dependent diabetes mellitus and non-insulin-dependent diabetes mellitus patients. The relationship to serum C-peptide immunoreactivity and free insulin. Ann Nutr Metab 1990; 34:333-342.
5. Lutz L, Barrett EJ, Groop LC, Ferrannini E, DeFronzo RA. Metabolic effects of low-dose insulin therapy on glucose metabolism in diabetic ketoacidosis. Diabetes 1988; 37:1470-1477.
6. Hale PJ, Crase J, Natrass M. Metabolic effects of bicarbonate in the treatment of diabetic ketoacidosis. Br Med J 1984; 289: 1035-1038.

**Distributed by:**  
Abbott Laboratories Ireland, Ltd.  
Abbott Diabetes Care  
4051 Kingswood Drive  
Citywest Business Campus  
Dublin 24  
Ireland  
1800 77 66 33

Abbott Diabetes Care Ltd.  
Abbott Diabetes Care  
Abbott House  
Yarnwall Business Park  
Yarnwall Road  
Maidenhead  
Berkshire  
SL6 4UD  
UK  
0500 467 466

FreeStyle and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions.  
This product(s) and/or its manufacture and/or use are protected by one or more of the following patents: US5,509,410; US5,628,890; US5,727,548; US6,129,823; US6,736,957; US6,764,581; US6,939,450; US6,377,894; US6,600,997; US6,773,671; US5,682,884; US6,591,125; US7,058,437; US7,504,019; US7,740,581; US7,905,999; US7,922,883; US7,998,337; US8,118,993; US8,182,671; US8,221,612; US8,211,280; US8,241,485; US8,241,486; US8,372,261; EP1,009,850B1; EP1,119,637B8; EP1,023,455B1; EP1135679B1; EP1801229B1; CA230248C; CA2346415C; CA2351796C; CA2353670C; CA2305800C. Additional patents may be issued and/or pending.

